


SIEMENS

Siemens Medical Solutions USA, Inc.
Mountain View, CA 94043, USA



Model: C5-2



For use with: ACUSON Freestyle



(1P) Model No.: 11001403 IPX8




(S) Serial No.: XXXXXX



FCC ID: XSB2300A
IC: 267FA-2300A


Contains FCC ID:
ED9LMX9838

Contains IC:
1520A-LMX9838




UWB: 屋外使用禁止

R 209 J00101



Siemens Medical Solutions USA, Inc.
Plymouth Meeting, PA 19462, USA



IVK

C5-2

SIEMENS

Siemens Medical Solutions USA, Inc.
Mountain View, CA 94043, USA



Model: L8-3



For use with: ACUSON Freestyle



(1P) Model No.: 11001401 IPX8



(S) Serial No.: XXXXXX




FCC ID: XSB2300A
IC: 267FA-2300A

Contains FCC ID:
ED9LMX9838

Contains IC:
1520A-LMX9838



R 209-J00099



UWB: 屋外使用禁止



Siemens Medical Solutions USA, Inc.
Plymouth Meeting, PA 19462, USA



IVK

L8-3

SIEMENS


Siemens Medical Solutions USA, Inc.
Mountain View, CA 94043, USA



CE
0123

Model: L13-5


For use with: ACUSON Freestyle



(1P) Model No.: 11001402 IPX8



(S) Serial No.: XXXXXX




07/2014

FCC ID: XSB2300A
IC: 267FA-2300A


Contains FCC ID:
ED9LMX9838

Contains IC:
1520A-LMX9838




R 209-J00100

UWB: 屋外使用禁止



Siemens Medical Solutions USA, Inc.
Plymouth Meeting, PA 19462, USA



IVK

L13-5