





FCC label

Bottom-case:ABS
Battery-cover:ABS
Shiny Surface
White

Nameplate



①34.0

②15.0

⚠ WARNING: Read and understand all instructions and warnings before use.
This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:
(1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operations.

4-R1.0③

AUTOMATIC BLOOD PRESSURE MONITOR
MODEL: KD-557BT Battery: 4×1.5V — SIZE AA Made in China
FCC ID: ZRY557BT
CE 0197    
ANDON HEALTH CO.,LTD.
NO.3 Jinping Street,YaAn Road,
Nankai District,Tianjin 300190,China SN